



PHC2001-8 Kombinerad pH-elektrod, röd stav, BNC

Produktnr: E16M313
SEK Pris: Kontakta oss
Tillgänglig

Long-term accuracy and reproducibility with Red Rod technology

General-purpose pH electrode with Red Rod reference system. Measuring range 0 to 12 pH, -10 to 100°C. Liquid junction: porous pin. With BNC plug.

The symmetry of Red Rod electrodes means that the iso pH is the same as the zero pH which provides highly reproducible results even when temperature fluctuations occur.

Reproducible measuring conditions ensured by a special bayonet head

Wide temperature range from -10 to 100°C

Fast and accurate response

Stable signal

Trouble-free liquid junction

Specifikationer

Elektrodtype:	PHC2001
Elektrolyt:	Refillable: Saturated KCl
Garanti :	12 månader
Givar koppling:	Porous pin
Givartyp:	Laboratory combination pH electrode with Red-Rod reference system
Kabellängd:	1 m
Längd:	103 mm
Material:	Glass
Mätområde:	0 - 12 pH
Specialfunktioner:	Red-Rod reference system
Temperaturområde :	-10 - 100 °C
Termistor:	None (external required)

Obligatoriska tillbehör

- Påfyllningslösning, referens, mättad KCl (KCl.L), 100 mL (Item S21M002)
- Påfyllningslösning, referens, KCl-kristaller (KCl.C), 15 g (Item S21M001)