



## pHC2701-8 Kombinerad pH-elektrod, högt KCl-flöde, BNC

<b>Produktnr:</b>	<b>E16M323</b>
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

### Long-term accuracy and reproducibility with Red Rod technology

pH electrode for high KCl flow with Red Rod reference system. Measuring range 0 to 12 pH, -10 to 100 °C. Liquid junction: annular ring. With BNC plug.

The symmetry of Red Rod electrodes means that the iso pH is the same as the zero pH which provides highly reproducible results even when temperature fluctuations occur.

### Reproducible measuring conditions ensured by a special bayonet head

Wide temperature range from -10 to 100°C

Fast and accurate response

Stable signal

Trouble-free liquid junction

---

## Specifikationer

Elektrodtyp:	pHC2701-8
Garanti :	6 månader