



## PHC2005-8 Robust kombinerad pH-elektrod, r d stav, BNC

**Produktnr:**

**E16M500**

SEK Pris:

Kontakta oss

Tillg nglig

### Long-term accuracy and reproducibility with Red Rod technology

General-purpose pH electrode with Red Rod reference system. Measuring range 0 to 12 pH, -10 to 100 C. Liquid junction: porous pin. With BNC plug.

The symmetry of Red Rod electrodes means that the iso pH is the same as the zero pH which provides highly reproducible results even when temperature fluctuations occur.

### Reproducible measuring conditions ensured by a special bayonet head

**Wide temperature range from -10 to 100 C**

**Fast and accurate response**

**Stable signal**

**Trouble-free liquid junction**

---

## Specifikationer

Elektrodtyp:

PHC2005

Garanti :

12 m nader