



Vatten, avjoniserat och avmineraliserat, 500 mL

Produktnr: 27249
SEK Pris: Kontakta oss
Tillgänglig

Demineralised Type 2 Water

Specific Conductance: $\leq 1 \mu\text{S/cm}$,

Silica (SiO_2): $\leq 30 \mu\text{g/L}$,

Chloride (Cl^-): $\leq 20 \mu\text{g/L}$,

Hardness: Negative,

Consumption of KMnO_4 : Negative

Specifikationer

Beskrivning:	Water, Deionized
Plattform:	Lösning
Volym:	500 mL