



Indikatorlösning med klorofosfonazo

Produktnr:

2589549

Farligt

SEK Pris:

Kontakta oss

Skickas inom 2 veckor



För bestämning av ultralåga hårdhetsvärden, 4-1000 µg/L CaCO₃.

Convenience

Rapid Liquid methods utilize liquid reagents, bottletop dispensers and pour-thru cells to provide simplified testing.

Performance

Since the pour-thru cell provides a stationary optical system, Rapid Liquid methods typically provide better precision than methods that require movement of sample cells.

Efficiency

Rapid Liquid methods are ideal for high-volume testing.

Specifikationer

Antal tester :	250
Mätområde:	4 - 1000 µg/L CaCO ₃
Metodnamn:	Chlorophosfonazo
Metodnummer:	8374
Volym/förpackningsstorlek:	500 ml