



REF201 Universell referenselektrod, 7,5 mm, röd stav

Produktnr:	E21M009
SEK Pris:	Kontakta oss
Tillgänglig	

General-purpose Red Rod reference electrode with banana plug. Liquid junction: porous pin, salt bridge solution: sat. KCl. Temperature range -10 to 100°C.

Red Rod electrodes feature:

- Reproducible measuring conditions
- Fast and accurate response
- Wide temperature range
- Stable signal
- Trouble-free liquid junction means less risk of clogging

Specifikationer

Garanti :	12 månader
Referenselektrod:	REF201