



REF251 Universal referenselektrod, 12 mm, röd stav, dubbel anslutning

Produktnr:

E21M001

SEK Pris:

Kontakta oss

Skickas inom 2 veckor

General-purpose Red Rod reference electrode with banana plug. Liquid junction: porous pin/porous pin, salt bridge solution: sat. KCl. Temperature range -10 to 100°C.

Red Rod electrodes feature:

- Reproducible measuring conditions
- Fast and accurate response
- Wide temperature range
- Stable signal
- Trouble-free liquid junction means less risk of clogging

Specifikationer

Garanti : 12 månader

Referenselektrod: REF251