



Reagenssats, hårdhet (total), digital titrator, fullt mätområde

Produktnr:

2272000

Farligt

SEK Pris:

Kontakta oss

Skickas inom 2 veckor



För bestämning av totalhårdhet, 10-160 mg/l CaCO₃ med EDTA-metoden.

For Total Hardness determination by EDTA Digital Titrator titration with Calmagite indicator. Range: 10-160 and 100-4000 mg/L as CaCO₃. Set contains ManVer® 2 Indicator, Hardness 1 Buffer Solution and one cartridge each of 0.0800 M and 0.800 M EDTA.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikationer

Mätområde:	10 - 160 mg/L & 100-4000 mg/L CaCO ₃
Metodnamn:	EDTA Titration
Parameter:	Hardness (Total)
Specialfunktioner:	Reagent set, Digital Titrator