



Filter paper, glass fiber, extra thick, pore size 1 µm, diameter 47 mm, 100/pk

Produktnr:	2842100
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

For filtration of heavily contaminated samples. High particulate-holding capacity makes discs efficient depth filters and allows for filtration of large volumes of solutions. Pall Life Sciences.

Specifikationer

Diameter:	47 mm mm
Material:	Glass Microfiber
Porstorlek:	1 µm
Volym/förpackningsstorlek:	100