



RADIOMETER M241Pt ORP-metallelektrod för allmänt bruk

Produktnr: E31M001
SEK Pris: Kontakta oss
Inget leveransdatum anges

The Radiometer M241Pt Metal Electrode has a platinum plate sensor and is recommended for general purpose measurements. An external reference electrode is required. The M241Pt has a 1 meter fixed cable with banana plug.

The right electrode for your application.

To meet your analysis needs, the Radiometer Analytical electrode range features single, double, and combined metal electrodes with platinum, silver, gold, antimony, mercury and glassy carbon sensing elements. This variety of designs can be used for redox measurements, redox titrations, potentiometric techniques, and imposed current potentiometric titration.

Specifikationer

Elektrodtype:	M241Pt
Garanti :	6 månader
Vad ingår i leveransen:	Inkluderar: Radiometer Analytical M241Pt metallelektrod med intyg om överensstämmelse.

Vad finns i förpackningen

Inkluderar: Radiometer Analytical M241Pt metallelektrod med intyg om överensstämmelse.