



HQ30D LBOD kit (cont. HQ30D, LBOD, probe, 6x BOD bottles, meter stand)

Produktnr:

HQ30D.99.701000

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

Rugged portable meters for use in the field and plant.

The digital meter/electrode system combines reliability, flexibility and ease of use. Interchangeable Intellical electrodes are recognised automatically and store all relevant data. Rugged, nearly indestructible outdoor versions with multiple cable lengths enable measurements even at previously inaccessible locations.

Replaced by NEW HQ Series.

Specifikationer

Dataexport:	Download via USB connection to PC or flash memory device. Automatically transfer entire data log or as readings are taken. Power supply adapter required.
Displaytyp:	240 x 160 pixel
DO-elektrodkalibrering:	0 - 1 DO calibration
DO-mätning:	0.00 - 20.0 mg/L
DO-upplösning:	0.01 mg/L eller 0.1 %
Garanti :	36 månader
Givare ingår:	None
Givartyp:	Intellical Digital SMART Standard Laboratory or Rugged Field
Ingångar:	1
Innehåll:	No
Intern resultatslagringskapacitet:	500 results
Kalibreringsintervall/varningar/påminnelse :	2 timmar till 7 dagar
Kapslingens Skyddsklass IP:	IP67
Kit?:	Ja
Konduktivitet, noggrannhet:	± 0.5 in (1µS/cm - 200mS/cm) range
Konduktivitet, upplösning:	5 digits with 2 digits after decimal point
Mätning av konduktivitet:	0.01 µS/cm to 200 mS/cm
Mätområde Redox:	-1500 - 1500 mV
Modell:	HQ30D – Multi/1 kanal
mV-mätning :	-1500 to 1500 mV

mV-upplösning :	0.1 mV
Needed Cable Length:	1
Omgivningsförhållande: relativ luftfuktighet :	90 % (non condensing)
Omgivningsförhållande: temperatur :	0 - 60 °C (32 - 140 °F)
Parameter:	pH, mV, ISE, DO, Conductivity, TDS, Salinity, Resistivity, ORP, Temperature
Parameters:	BOD
	NA
	NA
PC data transfer software :	HQD Series Meter Data Transfer Utility
pH Accuracy:	0,002 pH
pH Electrode calibration:	4 point calibration
pH Measurement Range:	0 - 14 pH
pH-upplösning :	0.001 - 0.1 pH
Resistivitmätning:	2,5 Ωcm - 49 MΩcm
Salinitetsmätning:	0 - 42 g/kg
Sensor A:	LBOD10101 (1 m)
Sensor B:	NA
Sensor C:	NA
Sensors:	Sensor A: LBOD10101
	Sensor B: NA
	Sensor C:
TDS-mätning:	0,0 - 50,0 mg/L
Temperaturmätning:	-10 - 110 °C
Temperaturnoggrannhet:	± 0,3 °C
Temperaturupplösning :	0,1 °C
Type of measurement:	QC Lab
Upplösning:	0,5
Vad ingår i leveransen:	Meter package includes HQ30D Portable Meter; 4 AA batteries; quick-start guide and user manual; HQD Meter Stand; Borosilicate glass BOD Bottles & ground glass stoppers, 300 mL, 6/ pk.
	Intellical LBOD101 Luminescent Dissolved Oxygen (LDO) sensor for BOD measurement, DO Cap with its iButton (already mounted), Test certificate, and Basic User Manual.
Vikt :	0,323 kg
Visning av kalibreringskurvor :	Kalibreringsöversikt

Vad finns i förpackningen

Meter package includes HQ30D Portable Meter; 4 AA batteries; quick-start guide and user manual; HQD Meter Stand; Borosilicate glass BOD Bottles & ground glass stoppers, 300 mL, 6/pk. Intellical LBOD101 Luminescent Dissolved Oxygen (LDO) sensor for BOD measurement, DO Cap with its iButton (already mounted), Test certificate, and Basic User Manual.