



## HQ40D Multimeter, Intellical BOD Electrode, 1 m

Produktnr:

HQ40D.99.701000

### UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

### Rugged portable meters for use in the field and plant.

The ultimate meter/probe system offering flexibility and ease of use with interchangeable probes. Ideal for facilities with multiple users and testing needs.

All HQD meters feature an advanced yet simple user interface that does not require manuals or training to operate. HQD meters connect with smart Intellical probes and automatically recognize the testing parameter, calibration history, and method settings to minimise errors and setup time. The durable portable meter and optional rugged probes are designed to withstand years of use in the field.

Choose from 4 models, including the dual-input HQ40D and single-input HQ30D versions. These multi-parameter meters offer the ultimate in measurement versatility - eliminating the need for multiple meters.

Replaced by NEW HQ Series.

---

### Specifikationer

Dataexport:	Download via USB connection to PC or flash memory device. Automatically transfer entire data log or as readings are taken. Power supply adapter required.
Garanti :	36 månader
Givare ingår:	None
Givartyp:	Intellical Digital SMART Standard Laboratory or Rugged Field
Ingångar:	2
Innehåll:	No
Kit?:	Ja
Konduktivitet, upplösning:	0.01 $\mu$ S/cm with 2 digits
Modell:	HQ40D – Multi/2 kanaler
Needed Cable Length:	1
Omgivningsförhållande: temperatur :	0 - 60 °C (32 - 140 °F)
Parameter:	pH, mV, ISE, DO, Conductivity, TDS, Salinity, Resistivity, ORP, Temp
Parameters:	BOD
	NA

	NA
PC data transfer software :	HQD Series Meter Data Transfer Utility
Sensor A:	LBOD10101 (1 m)
Sensor B:	NA
Sensor C:	NA
Sensors:	Sensor A: LBOD10101
	Sensor B: NA
	Sensor C:
Type of measurement:	QC Lab
Vad ingår i leveransen:	Meter package includes HQ40D Portable Meter; 4 AA batteries; Power adapter, USB/DC power adapter for data transfer; quick-start guide and user manual; HQD Meter Stand; Borosilicate glass BOD Bottles & ground glass stoppers, 300 mL, 6/pk.Intellical LBOD101 Luminescent Dissolved Oxygen (LDO) sensor for BOD measurement, DO Cap with its iButton (already mounted), Test certificate, and Basic User Manual.

---

## Vad finns i förpackningen

Meter package includes HQ40D Portable Meter; 4 AA batteries; Power adapter, USB/DC power adapter for data transfer; quick-start guide and user manual; HQD Meter Stand; Borosilicate glass BOD Bottles & ground glass stoppers, 300 mL, 6/pk.Intellical LBOD101 Luminescent Dissolved Oxygen (LDO) sensor for BOD measurement, DO Cap with its iButton (already mounted), Test certificate, and Basic User Manual.