



LTLCONPT100T 12 mm pH Electrode with Top68 Connector

Produktnr:	LZX537
SEK Pris:	Kontakta oss
Tillgänglig	

Low temperature glass and electrolyte, single reference system, electrode for ion-weak and/or cold water (10 - 100 µS/cm)

Specifikationer

Diameter:	12 mm
Elektrolyt:	KCl, saturated KCl/AgCl (crystal)
Garanti :	6 månader
Givargänga:	Pg 13,5
Kabelanslutning :	TOP68 connector
Längd:	140 mm
Mätområde:	-20 - 50 °C
Membran:	Ring diaphragm, porous PTFE
Neddopningsdjup:	105 mm
Reference system:	Ag/AgCl
Specialkommentarer:	For drinking water with low conductivity
Temperaturgivare:	Pt100
Tryckområde:	3.5 bar