



SOTA PT100T pH electrode, 12 mm, TOP68 connector

Produktnr: LZX583
SEK Pris: Kontakta oss
Inget leveransdatum anges

Standard glass and reference system with head connector TOP68

Specifikationer

Applikation:	General use
Elektrolyt:	KCl/AgCl + KNO ₃ gel, suitable for high temperatures
Garanti :	6 månader
Givargänga:	Pg 13,5
Kabelanslutning :	TOP68 connector
Kontakt:	TOP68
Längd:	140 mm
Material:	Sensor: Glass
Mätområde:	0 - 80 °C
Membran:	Ring diaphragm, porous PTFE, double
Neddopningsdjup:	105 mm
Referenselektrod:	None
Temperaturgivare:	PT100
Tryckområde:	0 - 3.5 bar