




Pole mount kit for LDO sensor, stainless steel

Produktnr:

LZX914.99.31100

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

Tank rim mounting with immersion pipe designed to mount onsite the luminiscent dissolved oxygen sensor.

Specifikationer

Diameter:	34 mm
Längd:	2300 mm immersion depth adjustable
Material:	Stainless steel 1.4571
Process connection:	CPVC clamp and stainless steel mounting bracket for the installation on the railing (max. diameter 2") or concrete wall
Specialkommentarer:	For the installation in tanks and open channels