



Moisture analyser IR-35, Denver, 35g capacity, 1mg (0.01%) readability

Produktnr: SMIR35

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

with metal rod radiator

Time Preset: 0.1-99 min, Data Interface: RS 232C, Characteristics: Resultdisplay in % Moisture, % Dry Matter, % ATRO, g Remainder, Determinationmode: Automatical Sample Heating: Infraredheating, via metal rod radiator, Results can be printed with opt. printer as a short journal or GLP-protocol in 6 different languages.