



pH Electrode for 8362 Ultrapure Water Sensor

Produktnr: 08362=A=2000

SEK Pris: Kontakta oss

Skickas inom 2 veckor

Part of analogue pH sensor 8362, specially designed for pure water pH measurement.

Specifikationer

Applikation:	Ultrapure water
Garanti :	6 månader
Givargänga:	Pg 13,5
Kabelanslutning :	S7 connector
Kontakt:	S7 Connector
Längd:	120 mm
Material:	Sensor: Glass
Mätområde:	0 - 14 pH
Referenselektrod:	None
Specialkommentarer:	pure water pH electrode
Temperaturgivare:	PT100