



Beaker, polypropylene, low form, 600 mL

Produktnr:	108052
SEK Pris:	Kontakta oss
Skickas inom 2 veckor	

Griffin low form. For general laboratory use. Non-toxic, noncontaminating, and biologically inert. Autoclavable. Combines no drip pouring, durability, and translucency. Approximate volume scales molded in for small sizes. Sizes 250 mL and larger feature silk-screened graduations and double spout for easier pouring. Nalgene.

Specifikationer

Diameter, yttre:	110.7 mm with spout
Enhet:	1 each
Instrument:	All except DR/4000, PC II
Kapacitet :	600 mL
Material:	Polypropylene, Low Form
Metodnummer:	8006