



## Funnel, analytical, 254 mL approx. Volume

<b>Produktnr:</b>	<b>108370</b>
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Standard long stem funnel meets the needs of a wide variety of applications. Mounded ribbing prevents "air lock" commonly occurring during funneling and filtrations. A steep, 58-degree cone slope makes for fast and easy flow of either liquids or small solids. Polypropylene construction is resistant to chemicals and makes this funnel fully autoclavable for sterilization.

---

### Specifikationer

Volym/förpackningsstorlek:	1 each
----------------------------	--------