



Test kit, nitrate-nitrite, model NI-12

Produktnr:

1408100

SEK Pris:

Kontakta oss

Inget leveransdatum anges

Farligt



Simple, low-cost color disc method is the most accurate visual test kit method.

Specifikationer

Antal tester :	100
Kapslingstyp:	D
Mätområde:	0 - 40 mg/L NO ₃ ⁻ -N
Metodnamn:	Color disc/Cadmium reduction, Diazotization
Modell:	NI-12
Parameter:	Nitrogen, Nitrate/Nitrite
Plattform:	Color Disc
Smallest increment:	1 0.01
Standardmätområde:	0 - 40 mg/L NO ₃ ⁻ -N; 0.01 - 0.40 mg/L NO ₂ ⁻ -N
Test-kit modell:	NI-12**