



Digital titratorpatron, natriumhydroxid, 3,636 N

Produktnr:

1438001

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



For carbon dioxide determination by Digital Titration

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikationer

Antal tester :	40 - 100
Mätområde:	100 - 1000 mg/L CO ₂
Metodnamn:	Phenolphthalein
Metodnummer:	8205
Parameter:	Carbon Dioxide
Plattform:	Digital Titrator
Titrant:	3.636 N NaOH