



m-FC broth, dehydrated, 100 g

Produktnr: 1462426

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking water samples, Hach has prepared media and convenient apparatus to make testing easy and fast.

Specifikationer

Media Type:	m-FC Broth
Volym/förpackningsstorlek:	100 g