



Laboratory Bottles, Polypropylene, Autoclavable, Wide Mouth, 500 mL, 12/pk

Produktnr:	1472459
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Polypropylene bottles suitable for sample collection to chemical storage. Matched polypropylene caps have a valve-seal design for a perfect, leakproof fit. A protruding ring under the cap holds shrink-wrapping securely in place. Autoclavable 121 °C (Do not autoclave with caps in place). Available in Narrow Mouth or Wide Mouth.

Specifikationer

Body Material:	Plastic
Kapacitet :	500 mL
Material:	PP
Package size UOM:	12pk
Volym/förpackningsstorlek:	12 /pk