



Phosphorus, orthophosphate (reactive) test kit, model PO-14

Produktnr:

147500

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Simple, low-cost color disc method is the most accurate visual test kit method

Specifikationer

Antal tester :	100 C**
Kapslingstyp:	D
Mätområde:	Single Range
Metodnamn:	Color disc/Stannous
Modell:	PO-14
Parameter:	Phosphorus
Plattform:	Color Disc
Smallest increment:	0.1
	1
Standardmätområde:	0 - 4.4 mg/L PO ₄
Test-kit modell:	PO-14