



Hardness, iron, and pH test kit, model HA-62B

Produktnr:

183702

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Model HA-62B

Includes color disc and color comparator for colorimetric comparison of iron and pH, as well as a drop-count titrator for hardness determination. Test kit produces a wide range of pH tests.

Combines frequently measured parameters in one rugged kit

Uses color disc and drop count titration methods

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikationer

Antal tester :	100
Applikation:	Water Quality
Modell:	HA-62B
Parameter:	Hardness, Iron, and pH (4-10)
Plattform:	Color Disc