



Cyanuric acid test kit, model CY-3

Produktnr: 185102

SEK Pris: Kontakta oss

Inget leveransdatum anges

A stabilizing agent for chlorine, cyanuric acid is used in swimming pool treatment. Cyanuric Acid Test Kit measures the turbidity formed in proportion to Cyanuric Acid concentration, after the reagent is added. Contains 100 tests.

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikationer

Antal tester :	100
Applikation:	Pool and Spa
Kapslingstyp:	D
Mätområde:	20 - 100 mg/L
Metodnamn:	Turbidimetry
Modell:	CY-3
Parameter:	Cyanuric Acid
Smallest increment:	20
Test-kit modell:	CY-3