



Tubing, tygon, R-3603, 4.8 mm (I.D.)

Produktnr:	2074043
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Chemical-resistant, Formula R-3603 Tygon tubing with flex strength greater than rubber. Unaffected by oxidizing agents and compatible with virtually any laboratory chemical. Bends readily to sharp radius curves and will maintain strong grip on glass or metal. Highly polished bore prevents interior build-up while facilitating cleaning.

Specifikationer

Diameter, yttre:	5/16 in
Formula (Text):	R-3603
I. D.:	3/16 in
Längd:	50 ft
Material:	Tygon
Slangväggstjocklek:	1/16 in