



Tubing, nalgene 8000, vacuum, 9.5 mm (I.D.)

Produktnr:	2074146
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Similar to Nalgene 8000 regular tubing, but with heavier walls. Imprinted in red for instant identification. Saves time and avoids errors which could invalidate experiments or injure laboratory personnel.

Specifikationer

Diameter, yttre:	7/8 in
Formula (Text):	R-3603
I. D.:	3/8 in
Längd:	50 ft
Material:	Tygon
Slangväggstjocklek:	1/4 in