



Bottle, BOD, black, 300 mL

Produktnr:	2086047
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

The standard for BOD analysis. Borosilicate glass bottle with flared lip and ground-glass stopper to prevent air entrapment and for water seal. Frosted marking area. Stopper included (60 and 300 mL sizes). Serialized bottles have silk-screened, sequential

Specifikationer

Body Material:	Glass
Kapacitet :	300 mL
Material:	Borosilicate Glass
Package size UOM:	1pk
Volym/förpackningsstorlek:	1 each