



Be Right™



Flask, Erlenmeyer, glass w/screw cap, 6/pk, 1000 mL

Produktnr: 2089683

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

Specifikationer

Enhet:	6/pk
Kapacitet :	1000 mL
Material:	Glass Screw Cap