



Be Right™



Flask, Erlenmeyer, polymethylpentene capacity 250 mL

Produktnr:	2089846
SEK Pris:	Kontakta oss
Tillgänglig	

Polymethylpentene flask with polypropylene screwcap features clarity, chemical resistance, strength, and autoclavability (Do not engage threads during cycle). Useful for titrations, preparation and storage of media.

Specifikationer

Enhet:	1 each
Kapacitet :	250 mL
Material:	Polymethylpentene
Typ:	Erlenmeyer