



Tube, culture, w/ screw cap, 20 x 150 mm, 10/pk

Produktnr:	2096230
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Cap has a firmly bonded rubber liner, providing a sure and highly inert sealing surface.

Specifikationer

Diameter:	20 mm
Dimension (D x L):	20 mm x 150 mm
Kapacitet :	Capacity to Neck 30 mL
Längd:	150 mm
Specialfunktioner:	Tube, 20 x 150 mm
Volym/förpackningsstorlek:	10 /pk