



Crucible, Gooch, 40 mL Capacity

Produktnr:	2098641
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Perforated bottom, glazed inside and out except for outside bottom surface. (Perforations are 1 mm in diameter). Filter for Gooch Crucible available and sold separately.

25 mL Gooch Crucible (# 2098640) is recommended for use in Standard Methods 2540 C and 2540 D to determine Total Suspended Solids.

Specifikationer

Kapacitet :	40 mL
Volym/förpackningsstorlek:	1 /pk