



Reagent, DPD total chlorine, test 'n tube vials, pk/50

Produktnr:

2105645

Farligt

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.



For total chlorine determination by the DPD Method. Test 'N Tube Reagent Vials contain pre-measured DPD reagent. Range: 0.09 - 5.00 mg/L Cl₂. Sample Size: 10 mL

Pack of 50 vials.

Test 'N Tube products provide safe, convenient testing

Capped 16-mm vials provide a self-contained package for mixing and measurement.

Easy ordering

All necessary reagents and vials are contained in the package.

Boric Acid free!

Specifikationer

Antal tester :	50
Instrument:	All except PC II
Mätområde:	0.09 - 5.00 mg/L Cl ₂
Metodnamn:	DPD
Metodnummer:	10101
Parameter:	Chlorine, Total
Plattform:	Powder Pillows