



Reagenspulverkuddar, kiseldioxid, citronsyra, 10 mL, pk/1000

Produktnr:

2106228

Farligt

SEK Pris:

Kontakta oss

Skickas inom 2 veckor



For determination of high range silica by the silicomolybdate method (Hach method 8185), and low range silica by the heteropoly blue method (Hach method 8186). Pack of 1000 Powder Pillows.

Specifikationer

Parameter:	Silica
Plattform:	Powder Pillows
Volym/förpackningsstorlek:	pkg 1000