



Reagenspulverkuddar, HR kiseldioxid, molybdat, 10 mL, pk/1000

Produktnr:

2107328

SEK Pris:

Kontakta oss

Tillgänglig

For determination of high range silica by the silicomolybdate method (Hach method 8185) Range:1 - 100 mg/L. Sample Size: 10 mL sample. Pack of 1000 Powder Pillows.

Specifikationer

Antal tester :	1000
Metodnamn:	Silicomolybdate
Plattform:	Powder Pillows
Volym/förpackningsstorlek:	1000 /pk