



Centrifuge Tube, Conical Bottom, 10/pk

Produktnr:	2115839
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Centrifuge tube with conical bottom. Polycarbonate with molded-in 0.1mL graduations. For low-speed centrifugation up to 6000 x g. Autoclavable. The 15 mL size is also used with Settlometer Kit (P/N 21157-00). Supplied in packages of 10, without closures. Autoclavable. Nalge.

Specifikationer

Kapacitet :	15 mL
Längd:	120 mm
Volym/förpackningsstorlek:	10 /pk