



Cylinder, graduated, polymethylpentene, 25 mL ± 0.3 mL, 0.5 mL divisions

Produktnr: 217240

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

This cylinder features transparent, nontoxic, noncontaminating, and biologically inert material with excellent chemical resistance. Meets ASTM E1272 Class B accuracy requirements, ISO Standard 6706, and suitable for food and beverage use (PMP resin meets requirements of CFR 21, Part 177.1520). Calibrated "To Contain/To Deliver" at 20°C (68°F). Shatterproof and thermally resistant to 135°C (275°F). Absence of meniscus for easy-to-read levels on front-to-back raised graduations.

Specifikationer

Certifieringar :	Class B
Kapacitet :	25 mL
Material:	Polymethylpentene
Mätmetod:	To Deliver/To Contain
Typ:	Single Metric
Volym/förpackningsstorlek:	1 each