



Dish, Petri, 15 x 100 mm

Produktnr:	2178900
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Widely used in microbiology, these standard-size polystyrene petri dishes are presterilized, stackable, and disposable. Dishes are packaged 24 per sleeve, with a case containing 25 sleeves (600 dishes total).

Standard-size polystyrene petri dishes

Presterilized, stackable, and disposable

600 dishes total

Specifikationer

Volym/förpackningsstorlek:	600
----------------------------	-----