



Test kit, formaldehyde, model FM-1

Produktnr:

2183100

SEK Pris:

Kontakta oss

Inget leveransdatum anges

Farligt



Formaldehyde solutions are used for maintenance of reverse osmosis systems and as disinfectants. The Model FM-1 measures higher concentrations by drop count titration. Contains 100 tests.

Utilises drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikationer

Antal tester :	100
Kapslingstyp:	D
Mätområde:	0,05 - 1% & 0,5 - 10%
Metodnamn:	Drop count titration/ Thymolphthalein
Parameter:	Formaldehyde - As CH_2O
Plattform:	Drop Count Titration
Smallest increment:	0.05% 0.5%
Test-kit modell:	FM-1