



Test kit, chromium, model CH-14

Produktnr:

222702

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Widely used as a corrosion inhibitor in cooling towers and engine cooling systems, chromium is also used in surface finishing, leather tanning, and other industries. The most widely used form, hexavalent chromium, is very toxic and strict limits are placed on discharges. Utilizing the drop count titration method to determine chromium (high range), this kit contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikationer

Antal tester :	100
Kapslingstyp:	D
Mätområde:	5 - 100 mg/L & 50 - 1000 mg/L Cr ⁶⁺
Metodnamn:	Drop count titration/ Thiosulfate
Parameter:	Chromium, high range - As Hexavalent Cr
Plattform:	Drop Count Titration
Smallest increment:	5
	50
Test-kit modell:	CH-14