



## Bag, sterile, WhirlPak, 207 mL, 25/pk

Produktnr:	2233125
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Complies with Standard Methods and USEPA for drinking water sample collection. Made from poly resins approved by USDA and FDA for food contact. Sterile (certified), disposable, transparent polyethylene bags with trusted "Whirl-Pak" closure. Suitable for liquid or solid sample collection. "Write-on" surface accepts waterproof ink. 24 oz bag perfect for incubating petri dishes in water bath applications.

---

### Specifikationer

Kapacitet :	207 mL
Volym/förpackningsstorlek:	25