



Test kit, chlorine (total), model CN-65

Produktnr:

225401

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



The most widely used disinfectant for drinking water, chlorine is also important for sanitizing swimming pools, cooling towers, other industrial equipment, and in the treatment of municipal wastewater. Its measurement and control are vital for both safety and economic reasons.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikationer

Antal tester :	100
Kapslingstyp:	D
Mätområde:	0,2 - 4,0 mg/L & 1 - 20 mg/L
Metodnamn:	Drop count titration/ Thiosulfate
Minsta ökningssteg:	1
Parameter:	Chlorine, mid and high range - As Total Cl ₂
Plattform:	Drop Count Titration
Smallest increment:	0.2
	1
Test-kit modell:	CN-65