



SA filter, cone, diameter 110 mm, soil analysis grade, 500/pk

Produktnr:	2263356
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Designed for sandy soils, SA filters are conically folded and have been optimized for use with the common extraction solutions used in soil nutrient analysis (Mehlich 1, Mehlich 3, Morgan, Modified Morgan, Bray, and Olsen). SA Filters provide the best combination of retention rating, filtration speed and purity requirements for soil nutrient testing (P, K, Ca, Mg). Medium filter retention.

Specifikationer

Diameter:	110 mm
Material:	Cellulose
Volym/förpackningsstorlek:	50 /pk