



Presence-absence (P-A) broth, dehydrated 500 g

Produktnr:

2280934

Farligt

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.



Used to detect pollution in treated or potable water; described in the 16th ed., Standard Methods.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of

your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make

wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking

water samples, Hach has prepared media and convenient apparatus to make

testing easy and fast.

Specifikationer

Media Type:	Presence-Absence (P-A) Broth
Volym/förpackningsstorlek:	500 g