



Be Right™



Bottle, Sampling, Polycarbonate, 1000 mL, 4/pk

Produktnr:	2324663
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Transparent polycarbonate bottle with cap features wide mouth (3.8 cm/1.48" diameter) opening for easy filling and square design for maximum space efficiency. Autoclavable with molded-in graduations (do not engage threads for autoclave cycle.)

Specifikationer

Body Material:	Polycarbonate
Kapacitet :	1000 mL
Volym/förpackningsstorlek:	4 /pk