



## Test kit, oxygen scavenger, color disc, model OS-1, 50 tests

**Produktnr:**

**2349500**

**Farligt**

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Proper maintenance and economy of operation of modern boilers and cooling towers require a careful balance between water quality parameters and treatment chemicals. Utilizing the color disc method, this kit contains 50 tests.

### Simple and fast color disc method

**Measures common oxygen scavengers: DEHA, hydroquinone, ISA, and MEKO**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

## Specifikationer

|                  |                                                                                             |
|------------------|---------------------------------------------------------------------------------------------|
| Antal tester :   | 50                                                                                          |
| Kapslingstyp:    | R                                                                                           |
| Mätområde:       | 0 - 300 µg/L DEHA (and other Oxygen Scavengers)                                             |
| Metodnamn:       | Color disc/Iron reduction                                                                   |
| Modell:          | OS-1                                                                                        |
| Parameter:       | Oxygen Scavenger - As DEHA,<br>hydroquinone,<br>erythorbic acid, and<br>methylethylketoxime |
| Plattform:       | Color Disc                                                                                  |
| Test-kit modell: | OS-1                                                                                        |