



Test kit, molybdate, model MO-LR

Produktnr:

2359300

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Molybdates are gaining over chromates for use as corrosion inhibitors due to lower toxicity. The high cost of these treatments warrant careful control of their concentration. Kit reagents react with Molybdate to form a yellow color, which is measured with a color disc.

Kit reagents react with Molybdate to form a yellow color, which is measured with an easy-to-use color disc method

Effectively monitors corrosion inhibitors

Specifikationer

Antal tester :	100
Kapslingstyp:	D
Mätområde:	Single Range
Metodnamn:	Color disc/ Ternary complex
Modell:	MO-LR
Parameter:	Molybdate
Plattform:	Color Disc
Smallest increment:	0.25
Standardmätområde:	0 - 3.0 mg/L Mo
Test-kit modell:	MO-LR