



Digital titratorpatron, kvicksilvernitratt, 2,570 N

Produktnr:

2393701

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



For determination of Salinity by Digital Titrator titration. 2.570 N. Digits x 0.1 = ppt Salinity in 2 mL sample.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikationer

Antal tester :	40 - 100
Instrument:	Digital Titrator
Mätområde:	0 - 100 ppt Salinity
Metodnamn:	Mercuric Nitrate
Parameter:	Sulfite
Plattform:	Digital Titrator
Titrant:	2.570 N Hg(NO ₃) ₂