



Test kit, nitrogen, ammonia, model NI-SA

Produktnr:

2428700

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Designed for testing ammonia in freshwater or seawater, this kit's color disc method offers low-cost, fast, reliable results in an easy-to-carry kit.

Designed for testing ammonia in freshwater or seawater

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifikationer

Antal tester :	100
Kapslingstyp:	D
Mätområde:	0 - 2.4 mg/L NH ₃ -N
Metodnamn:	Color disc/Salicylate
Parameter:	Nitrogen, Ammonia, mid range - For freshwater or seawater**
Plattform:	Color Disc
Smallest increment:	0.1
Test-kit modell:	NI - SA